# EV104

<b>GB</b>	<b>DIGITAL PERSONAL SCALE WITH A BMI INDEX</b>
<b>CZ</b>	<b>DIGITÁLNÍ OSOBNÍ VÁHA S BMI INDEXEM</b>
<b>SK</b>	<b>DIGITÁLNA OSOBNÁ VÁHA S BMI INDEXOM</b>
<b>PL</b>	<b>CYFROWA WAGA ŁAZIENKOWA Z POMIAREM BMI</b>
<b>HU</b>	<b>DIGITÁLIS SZEMÉLYMÉRLEG (BMI)</b>
<b>SI</b>	<b>DIGITALNA OSEBNA TEHTNICA Z INDEKSOM TELESNE MASE</b>
<b>RS HR BA</b>	<b>DIGITALNA LIČNA VAGA SA BMI INDEKSOM</b>
<b>DE</b>	<b>DIGITALE PERSONENWAAGE MIT BMI-INDEX</b>
<b>UA</b>	<b>ЦИФРОВИ ОСОБИСТІ ВАГИ З BMI ПОКАЗНИКОМ</b>
<b>RO</b>	<b>CÂNTAR DIGITAL DE PERSOANE CU INDICE BMI</b>
<b>LT</b>	<b>SKAITMENINĖS SVARSTYKLĖS</b>
<b>LV</b>	<b>DIGITĀLIE PERSONIESKIE SVARI</b>

# GB DIGITAL PERSONAL SCALE WITH A BMI INDEX EV104

Before using this product, study this manual and safety instructions carefully!

## Technical specifications:

Hardened safety glass: 6 mm  
 LCD display: 126 × 49 mm  
 Maximum weight: 150 kg  
 Minimum load: 2.5 kg  
 Measurement resolution: 0.1 kg/0.2 lb  
 Measurement units: kg/lb/st  
 Memory: max. 10 users  
 Height: 75–225 cm  
 Age: 10 to 85 years  
 Unit, body fat measurement range: 0.1%, 4.0%–50%  
 Unit, water content measurement range: 0.1%, 35%–85%  
 Unit, muscle mass measurement range: 0.1%, 10%–50%

Automatically turns off after 30 seconds  
 Overload/weak battery indication  
 Power supply: 2× 1.5V AAA (included)

## Basic rules of use:

- Digital personal scale PT-718 is designed solely for home use.
- Any repairs must be carried out by a professional service or by the seller!
- Do not use the scale near any appliances, which generate electromagnetic field. Put the scale far enough from such product so as to avoid interference.
- The result is only an indicative value. This product is not a medical appliance. If you want to reach optimum weight always consult your doctor before using any medicaments or diets.
- Do not put anything into the appliance's openings.
- The scale is designed for weighing of persons older than 10 years and younger than 85 years.
- Domestic animals and pets cannot be weighted on the scale.
- Protect the scale from falls and shocks. Do not jump on it.
- You or the scale can slip on wet surfaces. When standing on its smooth surface, never step on the edge and remember to maintain balance, when using the scale.
- Keep in mind the electronic impedance and avoid using this appliance in the following cases:
  - if you have a cardio stimulator or any other electronic implantation
  - if you have a cardiovascular disease
  - during pregnancy
  - if you go to a dialysis
  - if you are running a fever
  - if you suffer from osteoporosis
  - if you have an oedema or swelling
  - if you have been exercising for more than 5 hours
  - if your heartbeat is lower than 60 beats per minute.
- Use the scale only in compliance with the instructions from this manual.
- The producer is not responsible for damage caused by improper use of the appliance.
- This product is not to be used by persons (including children), whose physical, sensual or mental abilities or lack of experience and knowledge does not ensure safe use of the appliance, unless they are supervised or unless they have been instructed about the use of this appliance by a person in charge of their safety. Supervision over children is required to prohibit them from playing with the appliance.

## Error messages on the display

“8888” automatic test of functions, when the scale is turned on  
 “----” non-stable location of the scale, move it to a stable place  
 “Lo” Battery low – replace it.  
 Replace a depleted battery immediately to prevent damaging of the scale.  
 “Err” The body fat percentage is off the measuring range.  
 “EEEE” Scale overloaded, the maximum weight is 150 kg. Get off the scale to avoid its damaging.

## On screen icons

ICON	MEANING
kg (lb/st)	Unit of weight
USER	User memory
	Male
	Female
cm (ft/in)	Body height
Age	Age
	Non athletic build
	Average build
	Athletic build
	Fat content in body
	Water content in body
	Muscle content in body
	Bone content in body
%	Percentage of body fat / water / muscles / bones
	Indication of obesity: extremely thin – thin – normal weight – overweight – strongly overweight
KCAL	Daily calories intake unit
Increase	Weight increase since the last measurement
Decrease	Weight decrease since the last measurement

### Buttons:

SET – manually switch the scale on, set the input data

▲ – increase the value by 1 unit

▼ – decrease the value by 1 unit

### Battery replacement

1. Open the battery compartment at the bottom of the scale.
2. Remove the old batteries.
3. Insert new batteries. Maintain the proper polarity.
4. Close the cover.

### Basic measurement of weight

1. Put the scale on an even, hard and stable surface.
2. Carefully step on the scale and it will activate automatically.
3. Stand on the scale quietly, do not move. The measured weight will show up in a moment.
4. You can toggle the required weight units kg/lb/st by repeated pushing of the ▲ button.
5. After 30 seconds the scale will turn off automatically.

## MEASUREMENT OF BODY FAT, WATER AND MUSCLE MASS

### Measurement principle

Electric signals of a specific frequency pass through human body. The impedance of body fat is higher than the impedance of muscles or other tissues. The appliance generates an alternating electric current of a safe frequency and low voltage. Then the appliance analyses the biological resistance, which limits the passage of this current through human body (impedance of various tissues, fats and muscles differ) and from the different resistance there is calculated the percentage of fat in your body.

Fat % indicates the content of fat in human organism. Fat content in human body should vary within certain range. Too much body fat worsens metabolic processes and increases the risk of various diseases.

The displayed values are indicative only and cannot serve for medical purposes.

If you need a precise procedure how to achieve low weight, always contact a doctor.

### Measurement of weight with a display of body fat, water and muscle mass

Before weighing, you need to input basic data about each user.

You can store data about no more than 10 users.

Before weighing take off your shoes and socks, or the measurement will be incorrect.

### Push the SET button.

USER 0 (zero) sign will start blinking on the display.

Push SET again to store the information about the user number 0 or push ▲/▼ to change the user number.

Use ▲/▼ buttons to input information in the following order:

User number – sex – age – height – body built type.

To confirm each data, push SET.

You can change the stored data later, if you push SET.

After you enter all the data, the display will show 0.0 kg.

If you want to change the weight units, push ▲.

Step on the scale, stand still, do not move – after some time the following information will be displayed:

Current weight, fat content, water content, muscle mass, obesity degree icon →

Quantity of daily calories intake needed to maintain the current weight, bone content, muscle mass →

Decrease of weight (down arrow + Decrease) since the previous measurement or an increase of weight (up arrow + Increase).

This information will be shown one by one 3 times.

Then the scale will turn off.

### Maintenance and cleaning

1. To clean the scale use wet cloth. Make sure that no water gets inside. Keep the glass area dry. Do not use chemical cleaning aids with abrasive effect. Do not submerge the scale into water!
2. Keep the scale in horizontal position.
3. If the scale won't start, check if there is a battery and if it is not depleted.
4. If an error message is displayed or if the scale won't start, remove the battery for 5 seconds and reinsert.

Once the product's service life is over, do not throw it into unsorted communal waste. Use sorted waste collection points instead. Dispose of the used batteries in environmental friendly way and in compliance with your country's legislation. Proper liquidation of the product prevents negative impacts on human health and environment. Recycling of materials helps to save natural resources. You can get more information about the recycling of this product from your city council, your household waste disposal agency or the seller, where you purchased the product.

There is a certificate of conformity issued for this product.

13.8.2005

## OVERVIEW TABLE OF PHYSICAL CONDITION BASED ON THE FAT AND WATER CONTENT IN BODY

age	female		male		
	fat	water	fat	water	
<=30	4.0–20.5	66.00–54.7	4.0–15.5	66.00–58.1	<b>Low weight</b>
	20.6–25.0	54.6–51.6	15.6–20.00	58.00–55.00	<b>Normal weight</b>
	25.1–30.5	51.5–47.8	20.1–24.5	54.9–51.9	<b>Overweight</b>
	30.6–50.0	47.7–34.4	24.6–50.00	51.8–34.4	<b>Obesity</b>
>30	4.00–25.00	66.00–51.6	4.00–19.5	66.00–55.4	<b>Low weight</b>
	25.1–30.00	51.5–48.1	19.6–24.00	55.3–52.3	<b>Normal weight</b>
	30.1–35.00	48.00–44.7	24.1–28.5	52.2–49.2	<b>Overweight</b>
	35.1–50.00	44.6–34.4	28.6–50.00	49.1–34.4	<b>Obesity</b>

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